Instructions for use  
4BO/e™ RCM

mis®

### Resorbable Collagen Membrane 4BONE™ RCM

4BONE RCM is a bioresorbable collagen membrane, implantable material that is indicated for guided tissue regeneration procedure in periodontal defects to enhance regeneration of the periodontal apparatus.

#### DESCRIPTION

Resorbable collagen membrane 4BONE RCM is from highly purified collagen fibres type I and III derived from skin porcine origin. 4BONE RCM membrane is obtained by collagen cross linking using formaldehyde technology, and neutralized afterwards. 4BONE RCM membrane is fully resorbable eliminating the need for the second surgical procedure often required to remove a non resorbable membrane. Collagen is currently used for general and dental surgery as an absorbable haemostatic agent and absorbable wound dressing. 4BONE RCM is produced according to standardized, controlled manufacturing processes. 4BONE RCM resorbable membrane retards epithelial downgrowth during early phases of healing. Being semi occlusive, it allows essential nutrients to pass through. 4BONE RCM is very well incorporated to environment tissues and starts to resorb after 12 weeks. 4BONE RCM resorbable membrane has demonstrated excellent biocompatibility. 4BONE RCM is proposed in a sterile configuration. Opened, cannot be reused.

#### INDICATIONS

4BONE RCM is a bioabsorbable membrane. The 4BONE RCM is intended for use in periodontal/ dental surgery procedures as a material for placement in the area of periodontal defect, dental implant, bone defect or ridge reconstruction to aid in wound healing post surgery. Considering 4BONE RCM indications and resorption time, it is recommended to combine the membrane with a bone grafting product, for new bone healing by osteoconduction and osteoinduction e.g. 4BONE BCH.

#### PRECAUTIONS

The content of the double pouches is designed for single use only. Do not resterilize 4BONE RCM bone augmentation material. In case of membrane exposure during the healing phase the resorption time may be accelerated. If endosseous implants are involved, the membrane should be used only in combination with a stable implant and not in lieu of achieving primary implant stability. Absolute stability of the membrane is important for guided bone regeneration and is a vital condition for therapeutic success, and that the smallest movement on the tissue underneath is to be avoided.

#### INSTRUCTIONS

It is advisable that only clinicians trained in related treatment planning and in the technique of placing periodontal membranes should use 4BONE RCM absorbable collagen membrane. 4BONE RCM should be removed from the packaging using sterile gloves or instruments. Mucoperiosteal incision flaps are developed in the site to be treated. The incision should be as sulcular possible. Clinician should perform thorough debridement. As much as possible tissue must be preserved to allow for primary closure of the wound and correct positioning of the flaps. • Bone defect should be filled without excess with bone grafting material first, 4BONE BCH etc. • 4BONE RCM can be cut with sterile scissors to desired shape before being hydrated. 4BONE RCM can be placed either dry or hydrated. For better handling characteristics, the membrane can be hydrated in sterile water or saline solution for approximately not more than five minutes prior to the final placement. • 4BONE RCM should overlap the walls of the defect by at least 2mm to allow complete bone contact and to prevent gingival connective tissue invasion below the material. 4BONE RCM is applied over the defect without further treatment and held in place with moderate pressure. The period of time necessary to apply pressure will vary with the degree of bleeding. Adherence to the bone surface is achieved by gel formation of the collagen fibres with blood. • Complete penetration of the membrane by blood and exudates allows perfect adaptation and adhesion of the membrane to bony structures and makes the formation of a blood clot possible. • Salivary and other contamination to the material and surgical site should be minimized to avoid bacterial contamination. • Due to the high tensile strength of this membrane, fixation is possible. Fixation may be indicated to avoid membrane displacement due to loading or mobilization. • The mucoperiosteal flap is sutured over the collagen membrane (e.g. single sutures and deep mattress sutures). • The wound should be closed completely to avoid accelerated resorption due to membrane exposure.

### Special instructions for use in periodontology:

• A basic requirement for successful periodontal treatment includes eradicating the underlying bacterial

infection as well as adequate oral hygiene. Therefore, prior to surgical intervention, patients must receive a hygiene phase of treatment, consisting of oral hygiene instructions, scaling and root planning, and occlusal adjustment when indicated. A postoperative maintenance phase can help to ensure long-term therapeutic success. • In order to avoid the formation of a long junctional epithelium, 4BONE RCM must be adapted closely to the tooth (e.g. with additional fixation using suture material).

### LIMITATIONS FOR USE:

Due to the adherence to the bone tissue and the elasticity of 4BONE RCM, augmentation material is required to create and maintain space for bone formation. (e.g. 4BONE BCH) 4BONE RCM is therefore not indicated for single use in guided bone regeneration (GBR) without any space-making material. 4BONE RCM should not be placed where active infection exists. Before placement, the surgeon should be confident that any active or recent infection has been properly treated. The material has not been tested on pregnant women.

### POST OPERATIVE CARE

The patient should be monitored closely. • If the membrane becomes exposed, the dehiscence usually heals by itself within several weeks. Membrane removal is usually not necessary. However, to minimize bacterial contamination, rinsing with bactericidal solutions are recommended. • In the event that a membrane removal is necessary, the tissues adjacent to the membrane should be anesthetized with a local anesthetic. An incision should then be made immediately adjacent to the residual membrane. Following careful reflection of the surrounding tissue, the remaining portion of the membrane can be excised and the area curetted to remove any inflamed or infected tissue. To allow for undisturbed bone regeneration underneath the membrane surgical reentry should not be done before 4 to 6 months postoperatively.

### ADVERSE REACTIONS

4BONE RCM is high pure collagen membrane, adverse reactions after the use of 4BONE RCM have not been observed. Since 4BONE RCM is a collagen product, allergic reactions may not be totally excluded. Inflammation and adverse reactions like fever experienced with another microfibrillar collagen may be related. However, the following potential side effects may be noticed due to the surgical intervention; dehiscence, hematoma, pain, increased sensitivity and pain, redness, and inflammation.

### STERILISATION/CONSERVATION:

Do not use after expiry date. 4BONE RCM is sterilized by Beta irradiation, in double aluminum pouches. Do not resterilize. The consequences of a re-sterilisation and/or a re-use of the 4BONE RCM membrane have not been evaluated during the design process. So, the mechanical and resorption performances of 4BONE RCM membrane are not guaranteed in case of re-sterilisation and/or re-use. The membrane is sterile unless the package has been opened, damaged, or otherwise contaminated. Store at room temperature 15-25 °C (59°-77°F) and in a dry place.

### Key to codes used:

	Do not use if package is damaged		By prescription only
	Single use		Consult instructions for use
	Use by		Batch code
	Manufacture date (including sterilisation)		Sterilized using gamma irradiation
	Manufacturer		Catalogue number
			Federal law in USA and Canada restricts this device to sale by or on the order of a physician.

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MIS Implants Technologies Ltd.