Prime&Bond® NT

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Nano-Technology Etch&Rinse

Dental Adhesive System

CAUTION: This is a medical device. For dental use by dental professionals only.

USA: Rx only.

Content	Page
1 Product description	2
2 Safety notes	4
3 Step-by-step instructions	6
4 Hygiene	11
5 Lot number, expiration date and correspondence	12

1 Product description

Prime&Bond® NT Nano-Technology Etch&Rinse Dental adhesive is a self-priming dental adhesive system designed to bond resin based materials to enamel and dentin as well as to metals and ceramic.

Prime&Bond® NT adhesive combines primer and adhesive in a one-bottle material. The reduction of components and treatment steps simplifies use, maintaining superior bond strengths and protection against microleakage.

When mixed with Self Cure Activator, **Prime&Bond® NT** adhesive is designed to be used with Dentsply Sirona manufactured dual-cure/self-cure resin cements to bond all indirect restorations and Dentsply Sirona manufactured dual-cured composite restorations.

1.1 Intended purpose

Dental adhesive for tooth substrates and primer for restorative materials.

1.2 Clinical benefit

Restoration of tooth esthetics and function, as part of the restorative treatment.

1.3 Indications

Prime&Bond® NT adhesive (visible light-cured):

- Direct, light-cured composite and compomer restorations.
- Indirect restorations: light-cured, resin cemented veneers.
- Composite, ceramic and amalgam repairs.

Prime&Bond® NT adhesive mixed with Self Cure Activator (dual-cure):

- Direct, dual-cure or self-cure composite restorations and core build-ups.
- Indirect restorations: dual-cured and self-cured resin cemented inlays, onlays, crowns and bridge retainers.
- Dual-cured and self-cured resin cemented endodontic post cementation.

1.4 Contraindications

- Prime&Bond® NT adhesive is contraindicated for use in patients with a known hypersensitivity to methacrylate resins (e.g. Diethyleneglycol dimethacrylate (DGDMA)) or to any other component.
- Prime&Bond® NT adhesive is contraindicated for direct application to dental pulp tissue (direct pulp capping).

1.5 Delivery forms

Prime&Bond® NT adhesive is available in:

- Tri-layer bottle with flip-top-cap
- Unit Dose container 0.125 ml¹

1.6 Composition

Prime&Bond® NT adhesive:

- Acetone
- Urethane dimethacrylate (UDMA)
- Phosphoric acid modified acrylate resin (PENTA)
- Diethyleneglycol dimethacrylate (DGDMA)
- Elastomeric urethane dimethacrylate resin
- Functionalized amorphous silica
- Trimethacrylate
- Ethyl-4-dimethylaminobenzoate (photoaccelerator)
- Camphorquinone (photoinitiator)
- Butylated hydroxy toluene (BHT)
- Cetylamine hydrofluoride

1.7 Compatible restoratives/cements

- Prime&Bond® NT adhesive is chemically compatible with conventional (meth)acrylate-based visible light-cured composite restorative and cement materials including Dentsply Sirona visible light-cured composite and compomer restorative and cement materials.
- Prime&Bond® NT adhesive mixed with Self Cure Activator is chemically compatible with conventional (meth)acrylate-based dual and self-cured composite restorative and cement materials including Dentsply Sirona dual and self-cured composite restorative and cement materials (see complete Instructions for Use of selected restorative/cement material).

¹ Prime&Bond® NT adhesive in containers for single patient use is only available in selected countries.

2 Safety notes

Be aware of the following general safety notes and the special safety notes in other chapters of these Instructions for Use.

Safety alert symbol.

- This is the safety alert symbol. It is used to alert you to potential personal injury hazards.
- Obey all safety messages that follow this symbol to avoid possible injury.

Information symbol.

This information symbol is used to highlight important information for the correct use.

2.1 Warnings

The material contains polymerizable (meth)acrylate monomers, which may be irritating to skin, eyes and oral mucosa and may cause allergic contact dermatitis in susceptible persons.

This product contains Ethyl-4-dimethylaminobenzoate, which is classified as a CMR substance according to Regulation EC 1272/2008 on Classification, Labelling and Packaging of substances and mixtures. Based on a toxicological risk assessment, the presence of Ethyl-4-dimethylaminobenzoate does not pose a toxicological risk, if the product is used as intended.

- **Avoid eye contact** to prevent irritation and possible corneal damage. In case of contact with eyes, rinse with plenty of water and seek medical attention.
- **Avoid skin contact** to prevent irritation and possible allergic response. If contact with skin occurs, remove material with cotton and wash thoroughly with water and soap. In case of skin sensitisation or rash, discontinue use and seek medical attention.
- **Avoid contact with oral soft tissues/mucosa** to prevent inflammation. If accidental contact occurs, remove material from the tissues. Flush mucosa with plenty of water after the restoration is completed and expectorate/evacuate the water. If inflammation of mucosa persists, seek medical attention.
- For full coverage vital crown preparations, condition remaining enamel only. Etching of full coverage dentin surfaces is not recommended to minimize the possibility of post-operative sensitivity.

2.2 Precautions

This product is intended to be used only as specifically outlined in these Instructions for Use.

Any use of this product inconsistent with these Instructions for Use is at the discretion and sole responsibility of the dental practitioner.

- Safety and effectiveness in pregnant or breastfeeding women have not been established.
- Devices marked "single use" on the labeling are intended for single use only. Discard after use. Do not reuse in other patients in order to prevent cross-contamination.
- Use only suitable applicator tips (e.g. Applicator Tips by Dentsply Sirona) to ensure sufficient material is applied onto cavity surfaces.

- The bottles cannot be reprocessed. To prevent bottles from exposure to spatter or spray of body fluids or contaminated hands it is mandatory that the bottle is handled offside the dental unit with clean/disinfected gloves. During treatment clinicians with patient contact should not handle bottles.
- Contact with saliva and blood during application may cause failure of the restoration. Use of rubber dam or adequate isolation is recommended.
- Wear suitable protective eyewear, mask, clothing and gloves. Protective eyewear is recommended for patients.
- The Prime&Bond® NT adhesive bottles should be tightly closed immediately after use.
- Once Prime&Bond® NT adhesive is mixed with Self Cure Activator, apply immediately onto prepared surfaces.
- Proceed immediately once material is dispensed to avoid evaporation. Do not apply evaporated material.
- Use only in well ventilated areas. Avoid inhaling vapor.
- Flammable: Prime&Bond® NT adhesive and Self Cure Activator contain acetone. Keep away from sources of ignition. Take precautionary measures against static discharges.
- Interactions:
 - Eugenol- and hydrogen peroxide containing materials should not be used in conjunction with this product because they may interfere with hardening and cause softening of the polymeric components of the material.
 - Prime&Bond® NT adhesive is light-cured material. Protect from ambient light.
 - If mineral-impregnated (e.g., ferric compounds) retraction cords and/or hemostatic solutions are used in conjunction with adhesive procedures, marginal seal may be adversely affected, allowing microleakage, subsurface staining and/or restoration failure. If gingival retraction is necessary, use of plain, non-impregnated cord is recommended.
 - The use of a dual-cure adhesive system such as Prime&Bond® NT adhesive mixed with Self Cure Activator can shorten working time of a dual-cure resin cement system. This effect should be investigated in the laboratory prior to clinical use.
 - Various in vitro data exist regarding use of light-cured-only adhesives such as Prime&Bond® NT adhesive without Self Cure Activator in conjunction with most dual-cured or self-cured resin restoratives or cements. Chemical/Product incompatibility may adversely affect product efficacy, leading to premature restoration failure.

2.3 Adverse reactions

Product may irritate the eyes and skin.

- Eye contact: Irritation and possible corneal damage.
- Skin contact: Irritation or possible allergic response. Reddish rashes may be seen on the skin.
- Mucous membranes: inflammation, edema, sloughing (see [2.1 Warnings](#)).

Following side effects or adverse effects might be observed when using dental adhesives in general:

- Pulpitis resulting in post-operative pain or hypersensitivity.
- Degradation and failure of restoration due to marginal breakdown or debonding.
- Marginal discoloration.
- Secondary caries.

2.4 Storage conditions

Inadequate storage conditions may shorten the shelf life and may lead to malfunction of the product.

- Keep out of direct sunlight and store in a well ventilated place at temperatures between 10 °C-24 °C (50 °F-75 °F).
- Allow material to reach room temperature prior to use.
- Protect from moisture.
- Do not freeze.
- Do not use after expiration date.

3 Step-by-step instructions

3.1 General procedures

3.1.1 Cleaning: Preparation for all adhesive procedures

1. Use adequate isolation such as rubber dam.
2. Clean freshly instrumented enamel and dentin with water spray and then air dry. Clean uninstrumented enamel and dentin with a rubber cup and pumice or a cleaning paste such as Nupro® Prophylaxis Paste.
3. Wash thoroughly with water spray and blot or lightly air dry. Do not desiccate dentin.

3.1.2 Indirect restoration

1. Treat surface of restoration according to manufacturer's or dental laboratory's instructions, i.e., etching or mechanical roughening.
2. Apply Calibra® Silane Coupling Agent (available separately) to ceramic restorations designed to be silanated following manufacturer's instructions.

For zirconium, metal or other non-silanated restorations, follow adhesive application instructions below:

Light-cured cement, e.g., veneer cementation

1. Apply a single coat of Prime&Bond® NT adhesive to the internal bonding surface of the restoration. Immediately air dry for 5 seconds.
2. Light-cure Prime&Bond® NT adhesive for 10 seconds.²

Dual-/self-cured conventional resin cement, e.g., full coverage crown, endodontic post

1. Apply a single coat of Prime&Bond® NT adhesive mixed with Self Cure Activator to the internal bonding surface of the restoration.
2. Immediately air dry for 5 seconds.
3. Light-cure mixed adhesive/activator for 10 seconds.²

3.1.3 Repair of composite, ceramic and amalgam restorations with light-cured composites

1. Create mechanical retention where possible. Roughen the fractured surface for best results, sandblast the area to be repaired with an intraoral microetcher (50 µm alumina). Rubber dam is recommended with high speed evacuation.
2. Rinse with water for 10 seconds (microetched areas for 15-20 seconds), air dry.
3. For ceramic repair, etch with buffered hydrofluoric acid for intraoral use according to manufacturer's instructions and apply silane (e.g. Calibra® Silane Coupling Agent) to ceramic surfaces to be repaired following manufacturer's instructions.

² Curing light designed to cure materials containing camphorquinone (CQ) initiator. Peak of spectrum in the range of 440-480 nm and irradiance of 800-2000 mW/cm². Cure for at least 20 seconds if light output is between 550 and 800 mW/cm².

4. Proceed to pulp protection, tooth conditioning and adhesive application as outlined below.

3.2 Pulp protection

1. For direct and indirect pulp capping, cover the dentin close to the pulp (less than 1 mm) with a hard setting calcium hydroxide liner (Dycal® liner) leaving the rest of the cavity surface free for bonding with Prime&Bond® NT adhesive.

3.3 Tooth conditioning/dentin pretreatment

When used as a bonding agent for composite materials, it is recommended to follow the total etch/etch&rinse technique. When used as a bonding agent for compomer materials in non-stress bearing situations, acid etching is optional.

Prior to use, please see the complete Instructions for Use of the selected phosphoric acid conditioner.

3.3.1 Application of 34%-36% phosphoric acid conditioner

1. After application of rubber dam or other suitable isolation technique, apply the respective conditioner.
2. Apply the respective conditioner to the cavity surfaces starting at the enamel margins.
3. For best results, condition enamel for at least 15 seconds and dentin for 15 seconds or less.

3.3.2 Rinsing and drying

1. Remove conditioner with aspirator tube and/or vigorous water spray and rinse conditioned areas thoroughly for at least 15 seconds.
2. Remove rinsing water completely by blowing gently with an air syringe or by blot drying with a cotton pellet. Do not desiccate dentin.
3. Proceed immediately to application of Prime&Bond® NT adhesive.

Desiccation – To reduce the risk of inferior adhesion.

1. Do not desiccate the tooth substance – cavity surface should be adequate moist.
2. Avoid water pooling, leaving a moist, glistening surface. Do not rub the tooth surface when blot drying.

Contamination – To reduce the risk of restoration failure.

1. Once the surfaces have been properly treated, they must be kept uncontaminated.
2. If salivary contamination occurs, thoroughly clean with vigorous water spray, dry, and repeat conditioning procedure of enamel for 5 seconds only. Rinse and dry as described above.

3.4 Prime&Bond® NT adhesive application

Light-cured Procedures: Prime&Bond® NT adhesive.

- Direct, light-cured composite and compomer restorations.
- Indirect restorations: light-cured, resin cement veneers.
- Composite, ceramic and amalgam repairs.

Flip-top-cap Bottle

1. Pull off originality seal.
2. To open, take the bottle in the hand, place a thumb in the thumb recess and push the cap up until it remains in the open position.
3. Hold the bottle upside-down in a vertical position and dispense adhesive directly onto a clean Applicator Tip, making sure that the bottle does not come in direct contact with the brush, or place 1-2 drops into a standard dappen dish or mixing well.
4. If necessary clean dropper of the bottle using a soft paper tissue.
5. Carefully close the bottle again by applying pressure to the lever from above. The cap snaps into place when pressed closed with a perceptible "click".

Unit Dose container for single patient use

1. Grasp container at each end, or insert into holder, placing thumb along center score.
2. Firmly apply pressure until container separates.
3. Holder may be placed on the tabletop or held between fingers for convenience.
4. Insert disposable applicator into opening to saturate Applicator Tip.

Evaporation of solvent – To reduce the risk of inferior adhesion and sealing properties.

1. Immediately apply adhesive to maintain correct solvent concentration.

3.4.1 Application and curing

1. Using the disposable microbrush applicator tip supplied, immediately apply generous amounts of Prime&Bond® NT adhesive to thoroughly wet all the tooth surfaces. These surfaces should remain fully wet for 20 seconds and may necessitate additional applications of adhesive. Re-wetting of the microbrush may be required in order to coat the preparation for the full 20 seconds.
2. Remove excess solvent by gently drying with clean, dry air from a dental syringe for at least 5 seconds. Recommended technique is to begin drying/evaporation with the syringe approximately 10-15 cm (4-6 inches) away from the surface, gradually bringing the source to within 1 cm (0.5 inches) of the surface. Surface should have a uniform glossy appearance. If not, repeat application and air dry. Surface should not show areas of excessive adhesive thickness or pooling. Repeat air drying/evaporation outlined above if necessary.

Insufficient evaporation of solvent – To reduce the risk of inadequate polymerization.

1. Follow strictly the instructions of the steps above.

3. Cure Prime&Bond® NT adhesive for 10 seconds using a curing light designed to cure materials containing Camphorquinone (CQ) initiator. Peak of spectrum in the range of 440-480 nm and irradiance of 800-2000 mW/cm². Cure for at least 20 seconds if light output is between 550 and 800 mW/cm². Proceed immediately to completion/ restorative material placement.

Insufficient curing – To reduce the risk of inadequate polymerization.

1. Check compatibility of curing light.
2. Check curing cycle.
3. Check curing output before each procedure.

3.4.2 Completion

1. Place light-cured restorative material or light-cured resin cement over the cured Prime&Bond® NT adhesive as per restorative material manufacturer's instructions.
2. Ensure that the cured adhesive surfaces remain uncontaminated before applying the restorative material or cement. If salivary contamination occurs, thoroughly clean with vigorous water spray, dry and reapply adhesive as described above without repeating the acid conditioning step.

3.5 Prime&Bond® NT adhesive mixed with Self Cure Activator

Dual-cured or self-cured Procedures: Prime&Bond® NT adhesive with Self Cure Activator mixture.

- Direct, dual-cure or self-cure composite restorations and core build-ups.
 - Indirect restorations; dual-cured and self-cured conventional resin cemented inlays, onlays, crown and bridge retainers.
 - Dual-cured and self-cured resin cemented endodontic post cementation.
1. Clean preparation and treat restoration following 3.1 General procedures. Place Dycal® liner if needed following section 3.2 Pulp protection.
 2. For direct core build-up, place pin, post or matrix as needed.

3.5.1 Tooth conditioning/dentin pre-treatment

1. For full coverage crowns on vital teeth it is recommended not to etch dentin to minimize the possibility of post-operative sensitivity. Proceed to application of Prime&Bond® NT Dual-Cure Universal Dental Adhesive System, step 3.5.2 Application of Prime&Bond® NT adhesive mixed with Self Cure Activator.
2. When used as a bonding agent for dual-cure composites, resin cemented inlays, onlays, and endodontic posts, it is recommended to follow the total etch technique described in section 3.3 Tooth conditioning/dentin pretreatment.

Endodontic post space.

Use paper points to ensure removal of residual moisture.

3. Once the surfaces have been properly treated, they must be kept uncontaminated. If salivary contamination occurs, repeat procedure beginning at step 3.1 General procedures.

Desiccation – To reduce the risk of inferior adhesion.

1. Do not desiccate the tooth substance – cavity surface should be adequate moist.
2. Avoid water pooling, leaving a moist, glistening surface. Do not rub the tooth surface when blot drying.

Contamination – To reduce the risk of restoration failure.

1. Once the surfaces have been properly treated, they must be kept uncontaminated.
2. If salivary contamination occurs, thoroughly clean with vigorous water spray, dry, and repeat conditioning procedure of enamel for 5 seconds only. Rinse and dry as described above.

3.5.2 Application of Prime&Bond® NT adhesive mixed with Self Cure Activator

1. Place 1-2 drops of Prime&Bond® NT adhesive into a clean plastic mixing well. Close the bottle again promptly.
2. Place an equal number of drops of Self Cure Activator into the same mixing well. Replace cap promptly. Mix contents for 1-2 seconds with a clean, unused brush tip.
3. Using the disposable microbrush applicator tip supplied, immediately apply generous amounts of mixed Adhesive/Activator to thoroughly wet all the tooth surfaces. These surfaces should remain fully wet for 20 seconds and may necessitate additional applications of mixed adhesive/activator. Re-wetting of the microbrush may be required in order to coat the preparation for the full 20 seconds.

Endodontic post space.

Use of Endo-applicator tip is recommended. Alternatively, a paper point pre-wet with the adhesive mixture may aid in bringing the adhesive mixture down to the deepest portion of the preparation.

4. Remove excess solvent by gently drying with clean, dry air from a dental syringe for at least 5 seconds. Recommended technique is to begin drying/evaporation with the syringe approximately 10-15 cm (4-6 inches) away from the surface, gradually bringing the source to within 1 cm (0.5 inches) of the surface. Surface should have a uniform glossy appearance. If not, repeat application and air dry. Surface should not show areas of excessive adhesive thickness or pooling. Repeat air drying/evaporation outlined above if necessary.

Endodontic post space.

Use of clean, dry paper points may aid in thorough removal of solvent/ excess adhesive in post space.

Insufficient evaporation of solvent – To reduce the risk of inadequate polymerization.

1. Follow strictly the instructions of the steps above.

5. Light-cure mixed Adhesive/Activator for 10 seconds.³

Insufficient curing – To reduce the risk of inadequate polymerization.

1. Check compatibility of curing light.
2. Check curing cycle.
3. Check curing output before each procedure.

3.5.3 Completion

1. Dispense, mix and apply per manufacturer's directions:
 - Dual-cure or self-cure direct composite restorative/core build-up material.
 - Dual-cured resin cement for indirect restoration cementation.

³ Curing light designed to cure materials containing camphorquinone (CQ) initiator. Peak of spectrum in the range of 440-480 nm and irradiance of 800-2000 mW/cm². Cure for at least 20 seconds if light output is between 550 and 800 mW/cm².

4 Hygiene

Cross-contamination – To reduce the risk of infection.

- Do not reuse single use products. Dispose of in accordance with local regulations.
- Reprocess reusable products according to the instructions.

4.1 Bottles

Cross-contamination – To reduce the risk of infection.

1. The bottles cannot be reprocessed. Dispose of contaminated bottles in accordance with local regulations.
2. To prevent bottles from exposure to spatter or spray of body fluids or contaminated hands it is mandatory that the bottle is handled offside the dental unit with clean/disinfected gloves.
3. During treatment clinicians with patient contact should not handle the bottle.

Incidental contact of the bottle with water, soap or a water-based hospital-level disinfection solution will not damage bottle body. Do not allow any solution contact with contained material. Discard material that has been in contact with any fluid or non-sterile instrument.

Repeated liquid contact may damage label. Dry the bottle with a lint-free single-use cloth.

Vigorous wiping can destroy the label.

Wipe bottles gently.

4.2 Unit Dose Holder for Unit Dose container

For reprocessing instructions please refer to the Instructions for Use of the Unit Dose Holder, which is available on our webpage at <https://www.dentsplysirona.com/ifu>. If requested, we will send you a free printed copy of the Instructions for Use in the language you require within 7 days. For this purpose, use the order form provided on our webpage.

4.3 Disposal

Dispose of in accordance with local regulations.

5 Lot number (), expiration date () and correspondence

1. Do not use after expiration date.
ISO standard is used: “YYYY-MM” or “YYYY-MM-DD”
2. The following numbers should be quoted in all correspondence:
 - Reference number (REF)
 - Lot number (LOT)
 - Expiration date (EXP)
3. Any serious incident in relation to the product should be reported to the manufacturer and the competent authority according to local regulations.
4. Device Identification (Basic UDI-DI): ++D010ADH01MC

A summary of the safety and clinical performance (SSCP) for this product can be found (upon activation) at <https://ec.europa.eu/tools/eudamed> by searching using the Basic UDI-DI number listed above and at <https://www.dentsplysirona.com/ifu> using the reference number (REF).

For an explanation of the symbols used on the label, please refer to the Symbols Glossary at <https://www.dentsplysirona.com/ifu>.

© Dentsply Sirona 2025-08-11

[These Instructions for Use are based on Master Version 07]

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Dentsply DeTrey GmbH