Post length

Determine the post length and cut off the excess using only a diamond separating disc to avoid "fraying" the post and weakening it. Never use burs or clippers. We recommend cutting it outside the oral cavity. If you prefer maintaining the original length of the post during the reconstruction stages, it can always be shortened to the required length when core curing is completed. In this case it is also possible to use a normal bur.

Canal treatment

After using the endodontic instruments with the combined action of EDTA and Sodium hypochlorite for proper cleaning, and grinding with Universal Shaper burs, etching of the root canal dentin for 20 seconds is recommended.

Rinse and dry gently with paper cones and/or gentle air jets. Do not dehydrate.

Post insertion

Handle the post with suitable tools only and avoid contact with your hands. Should that be the case, clean the post with alcohol or a degreasing agent. Brush a dual enamel/dentine adhesive, following the manufacturer's instructions, onto both the canal walls and the post. Then use a composite cement, preferably self- light curing (dual) Bis-GMA-based as it is similar to the resin making up the post. The texture of the mixture must be fluid. Insert the cement into the canal using a syringe fitted with an endodontic nozzle and directly onto the post with a brush, in compliance with the manufacturer's instructions. Insert the post into the canal up to the chosen depth. If using a dual material, light cure the visible portion of the cement for 20 seconds, then wait for it to harden completely. If the material used in cementing is self-curing, wait until it has fully hardened (approx. 4 minutes). Do not use merely light-curing materials.

If using a dual cement – strongly recommended – insert the post as described above, but cure with a suitable light, moving the nozzle as close as possible to the fibre bundle of the post for at least 60 seconds. Then wait for it to harden completely (at least 5-6 minutes).

NOTE: Small posts may be used, at the doctor's discretion, under limited stress conditions (individually, e.g. in the incisal area, or in combination with other posts).

Core construction

The prosthetic core can be made using self or light-curing composites.

Packaging

See table.

For professional dental use only. Made in Italy.

ANATOMICAL POST - ANATOMICAL SUPER-LIGHT (TRANSLUCIDES) - ANATOMICAL POST S (STÉRILES)

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