DE**Gebrauchsanweisung****SUPRAMID** (PA 6-6/6)

Steriles, nicht-resorbierbares

Nahtmaterial**EN****Instructions for use****SUPRAMID** (PA 6-6/6)

Sterile, nonabsorbable

suture material**FR****Instructions d'utilisation****SUPRAMID** (PA 6-6/6)**Fil de suture**

stérile, non résorbable

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SUPRAMID

EN

Instructions for use, please read carefully!

The information given in this package leaflet is updated regularly. Please read instructions carefully prior to use.

POLYAMIDE 6-6/6

**Sterile, nonabsorbable
suture material**

DESCRIPTION

SUPRAMID is a pseudomonofilament nonabsorbable surgical suture material of polyamide 6/6, a polymer of hexamethylenediamine and adipic acid. It is coated with polyamide 6, a polymer of ϵ -caprolactam, in order to achieve an even surface. SUPRAMID is available either undyed (white) or dyed black (Carbon black, Colour Index 77266). SUPRAMID is available in various thicknesses and lengths and also as a needle-thread combination. SUPRAMID satisfies the requirements of the European Pharmacopoeia for sterile nonabsorbable suture material (0324).

INDICATIONS

SUPRAMID is used in general surgery and especially in ophthalmology. In addition, it is used in dental, implantological and cranio-maxillo-facial surgery. SUPRAMID is not used in neuro and cardiothoracic surgery.

MODE OF APPLICATION

The suture material to be used and its size should be chosen in accordance with the condition of the patient, the experience of the surgeon, the surgical technique, and the size of the wound. SUPRAMID may only be used by qualified medical personnel.

ACTION

SUPRAMID is used for coaptation of wound edges. It causes only a minimal inflammatory tissue reaction initially and is then slowly and progressively encapsulated by connective tissue. Although SUPRAMID is a nonabsorbable material, hydrolysis in vivo can lead to a progressive loss of tensile strength.

CONTRAINDICATIONS

As a progressive loss of tensile strength may occur after a prolonged period in vivo, SUPRAMID should not be used where tensile strength needs to be maintained over a prolonged period.

WARNINGS / PRECAUTIONS / INTERACTIONS

Care should be taken by known hypersensitivity against one of the

constituents of the nickel-chromium containing austenitic alloy of needles.

As with all suture material, care must be taken to ensure that the thread is not damaged during handling. In particular, it must not be kinked or crushed by surgical instruments such as needle holders.

When tightening the suture always pull on the thread between the needle and the puncture channel. Do not pull the thread too strongly or over sharp objects. When stretching the thread avoid friction with the rubber glove, as this can damage the thread.

In order not to damage the needle, always grasp it $\frac{1}{3}$ – $\frac{1}{2}$ of the distance from the end to the point.

Reshaping needles may cause them to lose strength and be less resistant to bending and breaking.

Users should exercise caution when handling surgical needles to avoid inadvertent needle stick injury. Discard needles in "Sharps" containers.

The risk of wound dehiscence varies with the location of the wound and the suture material used, therefore the user should be familiar with the surgical techniques in which SUPRAMID is to be used.

SUPRAMID must not be reused, there is a risk of infection for the patient and / or user and the functionality of the product can not be guaranteed.

ADVERSE REACTIONS

As with all foreign bodies, prolonged contact between SUPRAMID and saline solutions can lead to the formation of calculi (urinary and biliary tracts). Inflammatory tissue reactions and transient irritation in the area of the wound may occur initially. The foreign body effect can exacerbate preexisting tissue inflammation initially.

STERILITY

SUPRAMID is sterilised by ethylene oxide. Already opened or damaged packs must not be used! SUPRAMID is for single use only and should not be resterilized. Discard open, unused sutures! This is also applicable for Multiple-suture packs/ multipacks.

STORAGE CONDITIONS

Store at not more than 25°C protected from moisture and direct heat/ solar radiation. Do not use after the expiry date!

DESCRIPTION OF THE SYMBOLS USED ON THE PACKAGE

REF = Reference number **LOT** = Batch number

 = Use by year - month - day

 = Consult instructions for use

 = Do not reuse

 = Do not resterilise

 = Do not use if package is damaged

STERILE EO = Sterilised using ethylene oxide

 = Upper limit of temperature

 = Undyed, monofilament, nonabsorbable

 = Dyed, monofilament, nonabsorbable

 = Dyed, twisted, coated, nonabsorbable

PA = Polyamide  = HIBC code

 = Manufacturer

R_{only} = Prescription only (only for USA)

 = Package contents

 = Folding out the suture carrier to access the needle

CE 0197 = CE marking and identification number of the notified body. Product conforms to the essential requirements of the Council Directive 93/42/EEC concerning medical devices.

Needle code

1st letter: Indicates the curvature of the needle

F = 5/8-circle
H = 1/2-circle
D = 3/8-circle
V = 1/4-circle
K = Semi-curved
A = Asymptotic
G = Straight

2nd letter: Gives information on the type of needle and needle tip

R = Round-bodied needle
S = Cutting needle

3rd and 4th letters: Describe special features of the needle and needle tip

M = PREMIUM cut
(partly hand-honed)
N = Blunt point
T = Trocar point
P = Spatulated needle
S = Stronger diameter
X = Extra strong diameter
F = Thin PREMIUM cut "THIN LINE"
(partly hand-honed)
W = Flexible needle

Numbers indicate the straight (extended) length of the needle in millimetres

S (after numbers) = Stronger diameter
F (after numbers) = Extra thin diameter

H R X 30

Needle code

Example of needle size and needle curvature

= P = Spatula needle

= PM = Spatula needle PREMIUM-cut

= S = Reverse cutting needle

= SI = Inner cutting needle

= RT = Round-bodied cutting needle

= SM = Reverse cutting needle PREMIUM-cut

= SMI = Inner cutting needle PREMIUM-cut

= RTM = Round-bodied cutting needle PREMIUM-cut

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