Prior to use, carefully read the instructions for use.

EN

EPI-TEX™

FINISHING AND POLISHING STRIPS

For use only by a dental professional in the recommended indications.

INTENDED USE

For contouring, finishing and polishing of interproximal composite, compomer and glass ionomer restorative surfaces.

	Width	Length	Thickness	Grit	Use
Coarse (Blue)	5 mm	10 m	50 µm	#280	Contouring
Medium (Green)	5 mm	10 m	50 µm	#500	Finishing
Fine (Grey)	5 mm	10 m	50 µm	#800	Polishing
Extra Fine (Pink)	5 mm	10 m	50 µm	#1200	Polishing
Matrix strip (Transparent)	8 mm	10 m	50 µm	-	Contouring prior to light cure

1. Remove excess material with a suitable bur.
2. Choose a suitable grit.
3. Insert the strip between the interproximal contact point. Firmly grab both ends of the strip and draw back and forth.
4. Discard used strips.

CAUTION

- 1. EPI-TEX has sharp edges. Wear protective gloves and eyewear when handling the strips.
- 2. Take care not to injure the patient's soft tissues.
- 3. For infection control reasons, strips are for single use only.
- 4. Dispose of all waste according to local regulations.
- 5. Personal protective equipment (PPE) such as gloves, face masks and safety eyewear should always be worn.

PACKAGES

1. Starter Kit
1. refill (10 m) of each grit and Matrix Strip, 1 piece of dispenser
2. Refill pack available in all grits and Matrix strip, 1 refill (10 m)

STORAGE

Recommended for optimal performance, store away from high temperature or humidity.

Some products referenced in the present IFU may be classified as hazardous according to GHS. Always familiarize yourself with the safety data sheets available at: <http://www.gceurope.com> or for The Americas: <http://www.gcamerica.com>. They can also be obtained from your supplier.

Undesired effects- Reporting: If you become aware of any kind of undesired effect, reaction or other serious experience experienced by you with this product, including those not listed in this instruction for use, please report them directly through the relevant vigilance system, by following the proper authority of your country accessible through the following link: https://ec.europa.eu/growth/sectors/medical-devices/contacts_en, as well as to our internal vigilance system: vigilance@gce.dental In this way you will contribute to improve the safety of this product.

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