

mis[®]
PROSTHETIC OPTIONS
CONICAL CONNECTION

MIS superstructures are designed to meet all prosthetic needs and fit all established working protocols.

Abutments are available in a wide range of heights, gingival collar heights, angles and margin types. MIS offers solutions for all restoration types.





Table of Contents

Implant
Dimensions

P. 6 - 7

Prosthetic Options
Listed by Healing Caps

P. 8-16

Impression
Coping

P. 17-18

Temporary
Restoration

P. 19-21

Temporary
Cylinders

Screw-retained Restoration

P. 22-27

- CONNECT System
- Multi-Unit System
- Plastic Cylinders

Digital Restoration

P. 28-35

- EZ-Base System
- Ti-Base
- Digital Model Analogs
- Scan Posts
- Titanium Blanks

Cemented Restoration

P. 36-43

- Cementable Abutments
- CPK

Overdenture

P. 44-51


- OT-Equator Kits
- Ball Attachment System
- LOCKiT Attachment System
- Locator Attachment System

Prosthetic Screws


P. 52

Implant Dimensions: C1



Narrow
Platform



Standard
Platform






Wide
Platform







Implant Dimensions: V3

Narrow
Platform



Standard
Platform



NARROW PLATFORM

Prosthetic Options Listed by Healing Caps C1/V3



Each item chart displays available healing cap heights and related prosthetic options. **Bold item** numbers are shown in images above each column.

VN item numbers apply only to **V3 implants**.

1.

Healing caps

Impression coping - open tray



2.

Healing
capsImpression
coping -
open tray

or

Impression
coping -
closed trayCemented
options

or

Overdenture
options

3.

Healing
capsImpression
coping -
open tray



or

Impression
coping -
closed trayCemented
options

or

Screw
retained
options

or

Digital
restoration

STANDARD PLATFORM

Prosthetic Options Listed by Healing Caps C1/V3





Each item chart displays available healing cap heights and related prosthetic options. **Bold item** numbers are shown in images above each column.

VS item numbers apply only to **V3 implants**.

1.

Healing caps

Impression coping - open tray



CS-H0339	CS-I0395 CS-I0391
CS-H0439	
CS-H0539	CS-I0395
CS-H0639	
CS-H0839	

Gingival height 3, 4, 5, 6, 8mm

2.

Healing caps

Impression coping - open tray

or

Impression coping - closed tray

Cemented options

or

Screw retained options

or

Digital restoration

or

Overdenture options

3.



Healing caps

Impression coping - open tray

or

Impression coping - closed tray

Cemented options



WIDE PLATFORM

Prosthetic Options Listed by Healing Caps C1/V3

Each item chart displays available healing cap heights and related prosthetic options. **Bold item** numbers are shown in images above each column.

Healing
caps

Impression
coping -
open tray



or

Impression
coping -
closed tray

Cemented
options

Screw
retained
options

Overdentures



	CW-HS355	CW-IO500 CW-IO555 CW-IO501 CW-IO505	CW-IC551 CW-IC555 CW-IT300 CW-PF500	CW-AT001 CW-AN151 CW-SMAC1 CK-WPK61 CK-WPK81	CW-GP010 CW-GPC10 CW-W1480	CK-WOE1 CW-B1500
	CW-HS455			CW-AT002 CW-AT152 CW-AT252 CK-WPK62		
	CW-HS555		CW-IC551 CW-IC555 CW-PF500	CW-AT003 CW-AT153 CW-AT253 CW-MC103 CK-WPK63 CK-WPK64	CW-W3480	CK-WOE3 CW-B3500
Gingival height 3, 4, 5, 6mm	CW-HS655					


Healing
caps

Impression
coping -
open tray

or

Impression
coping -
closed tray

Digital
restoration



CW-HA263

CW-I0500
CW-I0555
CW-I0501
CW-I0505

CW-IC551
CW-IC555
CW-PF500
CW-IT300

CW-TB001
CW-TB002

CW-HA363

CW-TB015
CW-TBF15

CW-HA563

CW-IC551
CW-IC555
CW-PF500

CW-TB030
CW-TBF30

Gingival height 2, 3, 5mm

NARROW PLATFORM

Healing Caps V3/C1



Tightening
force:
20-25Ncm

Ø3.3

Height 2mm

3mm

4mm

5mm

6mm

8mm



CN-H0333



CN-H0433



CN-H0533



CN-H0633



CN-H0833

Ø4



CN-HS340



CN-HS440



CN-HS540



CN-HS640



CN-HS840

Ø4.8



CN-HA248



CN-HA348



CN-HA448



CN-HA548



CN-HA648



CN-HA848

STANDARD PLATFORM

Healing Caps V3/C1



Tightening
force:
20-25Ncm

Ø3.9

Height 2mm

3mm

4mm

5mm

6mm

8mm



CS-H0339



CS-H0439



CS-H0539




CS-H0639



CS-H0839

Ø4.8



CS-HS248



CS-HS348




CS-HS448



CS-HS548



CS-HS648



CS-HS848

Ø5.8




VS-HS358



VS-HS458



VS-HS558



VS-HS658



VS-HS858

WIDE PLATFORM

Healing Caps C1



Tightening force:
20-25Ncm

Height 2mm

3mm


4mm

5mm


6mm

8mm

Ø5.5



Ø6.3




Impression Coping

i Tightening force: 20-25Nm

Impression coping for open tray

Thin Impression coping for open tray



Impression Coping

i Tightening force:
20-25Ncm

Impression coping
for closed tray

Direct impression coping
for closed tray

Direct press fit for
closed tray

Implant analog

C1

V3

V3/C1

C1

NARROW

STANDARD

WIDE



CN-IC401




VN-IC401



CS-IC481



CW-IC551




CN-IC405




VN-IC405




CS-IC485




CW-IC555



CN-IT300




VN-IC400




CS-IT300




CW-IT300




CN-PF330



VN-PF340




CS-PF375




CW-PF500




CN-RSM10



VN-RSM10



CS-RSM10



CW-RSM10

Temporary Restoration Temporary Cylinders free-rotation

i Tightening
force:
20-25Ncm

Direct temporary
free-rotation cylinder

H1

Direct temporary
free-rotation cylinder


H2

Direct temporary
free-rotation cylinder

H3

V3/C1


NARROW



CN-TC013



CN-TC023



CN-TC033

V3/C1

STANDARD



CS-TC013



CS-TC023



CS-TC033



VS-TC013




VS-TC023




VS-TC033

C1

WIDE



CW-TC013



CW-TC023








CW-TC033

Temporary Restoration Temporary Cylinders anti-rotation

i Tightening
force:
20-25Ncm






Direct temporary
anti-rotation cylinder

H1







Direct temporary
anti-rotation cylinder

H2



Direct temporary
anti-rotation cylinder

H3



C1

V3

V3/C1

V3/C1

C1

NARROW

STANDARD


WIDE

Temporary Restoration Temporary Cylinders


i Tightening
force:
20-25Ncm

Direct temporary
plastic cylinder


H1




CN-TPA10



VN-TPA10



CS-TPA10



CW-TPA10

Direct temporary
plastic cylinder

H2




CN-TPH50



VN-TPH50



CS-TPH50



CW-TPH50

C1

V3

V3/C1

C1

NARROW

STANDARD

WIDE

Screw-Retained
Restoration
CONNECT
System

i Tightening
force:
30Ncm

H=1


H=2

H=3

H=4

V3/C1

NARROW




CN-C2040




CN-C3040

V3/C1

STANDARD



CS-C1540



CS-C2040




CS-C3040




CS-C4040

C1

WIDE



CW-C1540



CW-C2040



CW-C3040



CW-C4040


i

All CONNECT products
shown in this catalog are
sold subject to individual
regulatory certification only.

Screw-retained Restoration


CONNECT System Components

Healing cap




MM-H0540
MM-H1540
MM-H3040

Impression coping
for closed tray
Requires MT-IT100 key



MM-IC040

Impression coping
for open tray, free/
anti-rotation




MM-IO040
MM-IOI40

Abutment analog



MM-RSM41

Temporary cylinder,
free/anti-rotation



MM-TC041
MM-TCI41

Final esthetic abutment
free/anti-rotation



MM-CE046
MM-CEI46

Scan post
free/anti-rotation



MM-SP104
MM-SPF40

Prosthetic screw for
angled screw channel*
(Sold separately)



MM-SA160

CONNECT
prosthetic screw




MM-S0160

Long/short ratchet
insertion tool



MT-CLR21
MT-CSR21

Long/short motor
insertion tool



MT-CLM21
MT-CSM21

Digital model abalog



MM-MAN40

*For use with CONNECT final esthetic abutment
(MM-CE046, MM-CEI46) only and requires angled screw channel key (pg. 29)

Screw-retained
Restoration
Multi-Unit
System

i Tightening
force:
30Ncm

NARROW

C1

Height

1mm

2mm

3mm

4mm

5mm



CM-N1480



CM-N2480



CM-N3480



CM-N5480

STANDARD

V3/C1



CM-S1480



CM-S2480



CM-S3480



CM-S4480



CM-S5480

WIDE

C1



CM-W1480



CM-W3480

Screw-retained
Restoration
**Multi-Unit
System**
angulated
Abutments

i Tightening
force:
20Ncm

Height

1mm

2mm

NARROW


C1

STANDARD


V3/C1

WIDE


C1




CM-SO171




CM-SO301




CM-SO172




CM-SO302




CM-WO171



CM-WO301



CM-WO172



CM-WO302

Screw-retained Restoration


Multi-Unit System Components

Healing cap




MU-H4480

Analog



MU-RSM48

Impression coping,
open tray



MU-I0480

Impression coping,
closed tray



MU-IT480

Temporary cylinder




MU-TO480

Plastic cylinder




MU-CP048

Gold plastic cylinder



MU-BGC48
MU-MG480

CAD/CAM cementing
cap for multi unit



MU-CC480

Screw-retained Restoration Plastic Cylinders

i Tightening
force:
20-25Ncm

Direct anti-rotation
plastic cylinder

Direct free-rotation
plastic cylinder

C1

V3

V3/C1

C1

NARROW

STANDARD

WIDE



CN-GPC10




VN-GPC10



CS-GPC10



CW-GPC10



CN-GP010



CS-GP010



CW-GP020

Digital Restoration EZ-Base System

i Tightening
force:
20-25Ncm

EZ-Base
abutment

EZ-Base adjustable
abutment

C1

V3






V3/C1

C1

NARROW

STANDARD

WIDE



Digital Restoration

EZ-Base System Components




Tightening
force:
20-25Ncm

EZ-Base prosthetic screw




MD-E0220
MN-E0160

Long/short motor key for
angled screw channel



MT-ELR10
MT-ESR10

Long/short ratchet key for
angled screw channel



MT-ELM10
MT-ESM10




Unique screw and driver tip
design enable extreme angulation
and a wide path of insertion.


Digital Restoration
Ti-Base
anti-rotation

i Tightening
force:
30Ncm


H=0.5




CN-TB001



VN-TB001




CS-TB001




CW-TB001


H=1.5




CN-TB015



VN-TB015



CS-TB015



CW-TB015

H=3



CN-TB030



VN-TB030



CS-TB030





CW-TB030

Digital Restoration
Ti-Base
free-rotation

i Tightening
force:
30Ncm


H=0.5



H=1.5



H=3



Digital Restoration
Ti-Base
 Incisor esthetic
 anti-rotation

i Tightening
 force:
 30Ncm

H=0.5

H=1.5

H=3










C1

V3

V3/C1

NARROW

STANDARD



Digital Restoration Digital Model Analog



Tightening
force:
15Ncm

C1

V3

V3/C1

C1

V3/C1

V3/C1

NARROW

STANDARD

WIDE

MULTI-UNIT

CONNECT



CN-MAN10




VN-MAN10



CS-MAN10



CW-MAN10



MU-MAN10



MM-MAN40

Digital Dentistry Scan Posts



Tightening
force:
20Ncm

C1

V3

V3/C1

C1

V3/C1

V3/C1

NARROW

STANDARD

WIDE

MULTI-UNIT

CONNECT



CN-SP102



VN-SP102



CS-SP102



CW-SP102



MU-SP102



MM-SP104

Digital Restoration Titanium Blanks

i Tightening
force:
30Ncm

Ø9

Ø12

Amann Girrbach

Ø12

C1

V3

C1

C1

NARROW

STANDARD

WIDE



CN-IB091



VN-IB091



CS-IB091



CW-IB091



CN-IB121



VN-IB121



CS-IB121



CW-IB121



CN-AGT12



VN-AGT12



CS-AGT12



CW-AGT12



MN-MAB10



VN-MAB10



CS-MAB10







CW-MAB10

Cemented Restoration Cementable Abutments




i Tightening
force:
30Ncm

NARROW

0°





10°



15°

STANDARD

WIDE




Cemented Restoration Cementable Abutments

i Tightening
force:
30Ncm


NARROW

20°





25°

Cementing posts



STANDARD

WIDE



Cemented Restoration CPK - Transgingival Abutments

i Tightening
force:
30Ncm

C1

V3

V3/C1

C1

V3


NARROW

STANDARD


NARROW



CK-NPK41



VK-NPK41




CK-CPK41



CK-NPK61




VK-NPK61



CK-NPK42



VK-NPK42



CK-CPK42




CK-NPK62



VK-NPK62

C=4


C=6



CK-NPK43



VK-NPK43




CK-CPK43




CK-NPK63




VK-NPK63



CK-NPK44



VK-NPK44



CK-CPK44



CK-NPK64




VK-NPK64

STANDARD

WIDE

STANDARD




CK-CPK61

VK-CPK61




CK-WPK61




CK-CPK62

VK-CPK62



CK-WPK62




CK-CPK63

VK-CPK63



CK-WPK63




CK-CPK64

VK-CPK64




CK-WPK64

C=8




CK-CPK81

VK-CPK81




CK-CPK82

VK-CPK82



CK-CPK83

VK-CPK83



CK-CPK84

VK-CPK84

Cemented Restoration
CPK - Kit
Contents




Tightening
force:
30Ncm




Burnout plastic cap
anti-rotation




Burnout plastic cap
free-rotation



Plastic healing cap



Transfer
plastic cap



Analog

CPK abutments are
also sold in kits

CK-NPK41	●	Ø4, H=1, C=4
VK-NPK41	●	Ø4, H=1, C=4
CK-CPK41	●	Ø4.8, H=1, C=4
CK-NPK42	●	Ø4, H=2, C=4
VK-NPK42	●	Ø4, H=2, C=4
CK-CPK42	●	Ø4.8, H=2, C=4
CK-NPK43	●	Ø4, H=3, C=4
VK-NPK43	●	Ø4, H=3, C=4
CK-CPK43	●	Ø4.8, H=3, C=4
CK-NPK44	●	Ø4, H=4, C=4
VK-NPK44	●	Ø4, H=4, C=4
CK-CPK44	●	Ø4.8, H=4, C=4
CK-NPK61	●	Ø4, H=1, C=6
VK-NPK61	●	Ø4, H=1, C=6
CK-CPK61	●	Ø4.8, H=1, C=6
CK-NPK62	●	Ø4, H=2, C=6
VK-NPK62	●	Ø4, H=2, C=6
CK-CPK62	●	Ø4.8, H=2, C=6
CK-NPK63	●	Ø4, H=3, C=6
VK-NPK63	●	Ø4, H=3, C=6
CK-CPK63	●	Ø4.8, H=3, C=6
CK-NPK64	●	Ø4, H=4, C=6
VK-NPK64	●	Ø4, H=4, C=6
CK-CPK64	●	Ø4.8, H=4, C=6
CK-CPK81	●	Ø4.8, H=1, C=8
VK-CPK81	●	Ø5.5, H=1, C=8
CK-CPK82	●	Ø4.8, H=2, C=8
VK-CPK82	●	Ø5.5, H=2, C=8
CK-CPK83	●	Ø4.8, H=3, C=8
VK-CPK83	●	Ø5.5, H=3, C=8
CK-CPK84	●	Ø4.8, H=4, C=8
VK-CPK84	●	Ø5.5, H=4, C=8




Overdenture OT-Equators

i Tightening
force:
30Ncm


NARROW

V3/C1




STANDARD

V3/C1



WIDE

C1



Overdenture OT-Equator Components

i Tightening
force:
30Ncm

OT-Equator
handling tools



ET-IT002

OT-Equator impression
coping cap




OE-IC002

OT-Equator impression
transfer cap



OE-IC001

OT-Equator analog




OE-RSM60

Metal housing &
transfer caps



OE-MH001 OE-RCB01
OE-RCW01 OE-RCP01
OE-RCY01 OE-RCV01

OT-Equators
are sold in kits



Overdenture Ball Attachment System

i Tightening
force:
30Ncm

NARROW

V3/C1



CN-B1330



CN-B2330




CN-B3330



CN-B4330

STANDARD

V3/C1



CS-B1375



CS-B2375



CS-B3375




CS-B4375



CS-B5375

WIDE

C1



CW-B1500



CW-B3500

Overdenture Ball Attachment Components

i Tightening
force:
30Ncm

Analog




MB-RS375

Gold / Plastic /
Titanium ball cap



MB-GPS10 MB-TBS20
MB-SF200 MB-TBS10
MB-TTP10 MB-PR010

Metal housing &
Transfer caps



MB-DMH10
MB-PPC10
MB-PPR10

Plastic disc for ball
attachment



MB-DB235

Ball attachment
anchor key



MT-NB225

Ball attachment
anchor key



MT-RB225

Overdenture LOCKiT Attachment System

i Tightening
force:
30Ncm

NARROW

V3/C1



STANDARD

V3/C1



WIDE

C1



All LOCKiT products shown in this catalog are sold subject to individual regulatory certification only.

Overdenture LOCKiT System Components

i Tightening
force:
30Ncm

Abutment holder sleeve



MT-LOC11

Caps range




MM-BLA10 MM-CLE10
MM-BLU10 MM-RED10
MM-PIN10 MM-GRE10

Male processing package



MM-LOC02

Metal housing



MM-TIT10

Impression coping



MM-LOC1

LOCKiT analog



MM-LOA40
MM-LOA50

LOCKiT core tool




MT-LOC10

LOCKiT long ratchet
driver



MT-RSD24

LOCKiT short ratchet
driver



MT-RSD19

LOCKiT abutments
are also sold in kits


CK-NLC1 ● Ø3.75, H=1
CK-SLC1 ● Ø3.75, H=1
CK-WLC1 ● Ø4.5, H=1

CK-NLC2 ● Ø3.75, H=2
CK-SLC2 ● Ø3.75, H=2

CK-NLC3 ● Ø3.75, H=3
CK-SLC3 ● Ø3.75, H=3
CK-WLC3 ● Ø4.5, H=3

CK-NLC4 ● Ø3.75, H=4
CK-SLC4 ● Ø3.75, H=4

CK-NLC5 ● Ø3.75, H=5
CK-SLC5 ● Ø3.75, H=5
CK-WLC5 ● Ø4.5, H=5



Overdenture Locator Attachment System

i Tightening
force:
30Ncm

NARROW

V3/C1



STANDARD

V3/C1



Overdenture Locator System Components



Tightening
force:
30Ncm

Analogs



MM-LFA40
MM-LFA50

Metal housing &
transfer caps



MM-LMP02

Impression coping



MM-LIC10

Caps
range:



MM-LLR10



MM-LER10



MM-LELE1



MM-LRM10



MM-LEL10



MM-LPP10



MM-LBP10

Tools:



MM-LCT10



MM-RSD15



MM-LSD15



MM-RSD20



MM-LSD21

Prosthetic Screws



Tightening
force:
20/30Ncm

Mini prosthetic screw




MD-S0200

Direct prosthetic screw



MN-S0160
MD-S0220

Medium prosthetic screw



MD-S0222

Long prosthetic screw



MD-S0224

Direct gold prosthetic screw




MD-G0220

NOTES



NOTES





NOTES



