

mis[®]
PROSTHETIC OPTIONS
INTERNAL HEXAGON



Table of Contents

Implant
Dimensions

P. 8-9

Prosthetic Options
Listed by Healing Caps

P. 10-14

Analogs and
Impression Copings

P. 15-16

Temporary
Restoration

P. 17-19

Temporary
Cylinders

Screw-retained Restoration

P. 20-25

- Multi-Unit System
- Plastic Cylinders

Digital Restoration

P. 26-31

- EZ-Base System
- Ti-Base
- Digital Model Analogs
- Scan Posts
- Titanium Blanks

Cemented Restoration

P. 32-41

- Cementable Abutments
- CPK

Overdenture

P. 42-49

- OT-Equator Kits
- Ball Attachment System
- LOCKIT Attachment System
- Locator Attachment System

Prosthetic Screws


P. 50

Implant Dimensions: SEVEN





Narrow
Platform



Standard
Platform



Wide
Platform



Implant Dimensions: M4






Narrow
Platform



Standard
Platform




Wide
Platform



Prosthetic Option Listed by Healing Caps


The internal hex. prosthetic system includes the existing abutments. In addition, there is a line of consistent concave emergence profile abutments for the benefit of the entire rehabilitation process.






Each item chart displays available healing cap heights and related prosthetic options. **Bold item** numbers represent image shown above each column.

*Image is for illustration only. (2mm height closed tray impression coping presented. Not available in 3mm).

- Healing caps
- Impression coping - open tray
- or
- Impression coping - closed tray
- or
- Cemented options
- or
- CAD/CAM solutions





NARROWPLATFORM

Healing Caps



Tightening
force:
20-25Ncm

Healing cap

Ø4

HEIGHT 2mm



MH-N2330

3mm



MH-N3330

4mm



MH-N4330

5mm



MH-N6330

6mm

8mm




MH-N8330

Concave healing cap

Ø4



MH-N0440



MH-N0540



MH-N0640

Anatomic healing cap


Ø4.3



MH-52330




MH-53330



MH-54330



MH-55330



MH-56330


STANDARD PLATFORM Healing Caps

i Tightening
force:
20-25Ncm

HEIGHT 2mm 3mm 4mm 5mm 6mm 8mm


Healing cap

Ø4




Healing cap

Ø4.8




Concave healing cap

Ø4.8



Anatomic healing cap

Ø5.5



WIDE PLATFORM

Healing Caps



Tightening
force:
20-25Ncm

HEIGHT

2mm

3mm

4mm

5mm

6mm

8mm

Healing cap

Ø5



MH-W3500



MH-W4500



MH-W5500

Concave healing cap

Ø5.5



MH-W0455



MH-W0555



MH-W0655

Anatomic healing cap

Ø6.5



MH-W3630



MH-W4630



MH-W5630

Impression coping




Tightening force:
20-25Ncm

Impression coping for closed tray

NARROW

STANDARD


WIDE



MN-IT100



MD-IT100




MW-IT100


Direct impression coping for closed tray



MN-IT300




MD-IT300




MW-IT300

Direct impression coping for closed tray, concave profile



MN-ITC30



MD-ITC30



MW-ITC30

Temporary Restoration

Temporary Cylinders -




Tightening
force:
20-25Ncm

Narrow profile
temporary plastic
abutment


Temporary plastic
abutment

Direct temporary
plastic cylinder

NARROW




MN-TPA12



MN-TPA10

STANDARD




MD-TPA12




MD-TPA10

WIDE



MW-TPA12



MW-TPA10



MD-TPH50



MW-TPH50


Temporary Restoration
Temporary
Cylinders -
free-rotation

i Tightening
force:
20-25Ncm


Direct temporary
free-rotation cylinder

H1


NARROW



STANDARD




WIDE




Direct temporary
free-rotation cylinder,
concave profile

H2



Direct temporary
free-rotation cylinder,
concave profile

H3



Temporary Restoration


Temporary Cylinders - anti-rotation

i Tightening
force:
20-25Ncm


Direct temporary anti-
rotation cylinder

H1


NARROW



STANDARD



WIDE





Direct temporary
anti-rotation cylinder,
concave profile

H2



Direct temporary
anti-rotation cylinder,
concave profile

H3



Screw-retained
Restoration

Multi-Unit System



Tightening
force:
30Ncm



H=2

NARROW



MA-NC248

STANDARD



MA-SC248

WIDE

H=3



MA-NC348



MA-SC348



MA-WC348

H=4



MA-NC448



MA-SC448

H=5



MA-WC548

Screw-retained
Restoration

Multi-Unit System



Tightening
force:
30Ncm

H1

NARROW



MA-N1480

STANDARD



MA-S1480

WIDE



MA-W1480

H2



MA-N2480



MA-S2480


H3



MA-N3480



MA-S3480



MA-W3480

H4




MA-N4480



MA-S4480

H5



MA-N5480



MA-S5480



MA-W5480

Screw-retained
Restoration


Multi-Unit System



Tightening
force:
20Ncm

17°


STANDARD



MA-SO171



MA-SO172



MA-SO301




MA-SO302


WIDE




MA-WO172



MA-WO171



MA-WO301




MA-WO302

30°

Screw-retained Restoration


Multi-Unit Components

Healing cap




MU-H4480

Analog



MU-RSM48

Impression coping,
open tray



MU-I0480

Impression coping,closed tray



MU-IT480

Temporary cylinder




MU-TO480

Plastic cylinder



MU-CP048

Gold plastic cylinder



MU-BGC48
MU-MG480

CAD/CAM cementing
cap for multi unit



MU-CC480

Screw-retained
Restoration

Plastic Cylinders



Tightening
force:
20-25Ncm

Direct plastic
cylinder, anti-rotation

Direct plastic
cylinder free rotation


Wide profile direct plastic
cylinder, anti-rotation

Wide profile direct plastic
cylinder, free rotation

NARROW



MN-CPH13




MN-CP013

STANDARD



MD-CPH13



MD-CP013

WIDE




MW-CPH13



MW-CP013




MD-CPH50



MD-CP050



MW-CPH60



MW-CP060

Screw-retained Restoration

Plastic Cylinders

i Tightening
force:
20-25Ncm

Direct gold plastic
cylinder, anti-rotation

Direct gold plastic
cylinder, free rotation

Direct CoCr plastic
cylinder, anti-rotation

Direct CoCr plastic
cylinder, free-rotation

NARROW



MN-GPC10



MN-GP010



MN-CCI10



MN-CC010


STANDARD




MD-GPC10



MD-GP010



MD-CCI10



MD-CC010

WIDE



MW-GPC10



MW-GP010



MW-CCI10



MW-CC010


Digital Restoration EZ-Base System

i Tightening
force:
25Ncm


NARROW



STANDARD



WIDE



Digital Restoration EZ-Base Components

i Tightening
force:
25Ncm

EZ-Base prosthetic screw



*MD-E0220
*MN-E0160

* For use with ratchet/motor key for angled screw channel



MT-ELR10
MT-ESR10



MT-ELM10
MT-ESM10

Digital Restoration Ti-Base

i Tightening
force:
30Ncm

CONCAVE

Anti-rotation

Free-rotation


NARROW

STANDARD


WIDE




MN-TBC15



MD-TBC15



MW-TBC15




MN-TBC30



MD-TBC30




MW-TBC30



MN-FTB15



MD-FTB15



MW-FTB15



MN-FTB30



MD-FTB30






MW-FTB30


Digital Restoration Ti-Base

i Tightening
force:
30Ncm


Anti-rotation



Free-rotation



Incisor esthetic



Digital Restoration

Digital Model Analog



Tightening
force:
15Ncm

NARROW



MN-MAN10

STANDARD



MD-MAN10

WIDE



MW-MAN10

MULTI-UNIT



MU-MAN10

Digital Dentistry
Scan Posts



Tightening
force:
30Ncm

NARROW

STANDARD

WIDE

MULTI-UNIT



Digital Restoration

Titanium Blanks

i Tightening
force:
30Ncm

Amann Girrbach

NARROW

Ø9mm



MN-IB091

Ø12mm



MN-IB121

Ø12mm



MN-AGT12



MN-MAB10

STANDARD

Ø9mm



MD-IB091

Ø12mm



MD-IB121

Ø12mm



MD-AGT12



MD-MAB10

WIDE

Ø9mm



MW-IB091

Ø12mm



MW-IB121

Ø12mm



MW-AGT12



MW-MAB10


Cemented Restoration
Cementable
Abutments

i Tightening
force:
30Ncm



0°

NARROW


STANDARD












10°



15°




WIDE




Cemented Restoration
Cementable
Abutments

i Tightening
force:
30Ncm

20°




NARROW



25°



STANDARD





Cemented Restoration Cementing Posts

i Tightening
force:
30Ncm



NARROW



STANDARD



WIDE



Cemented Restoration


CPK - Transgingival Abutments

i Tightening
force:
30Ncm


NARROW

STANDARD


WIDE




MN-CPK41




MD-CPK41




MW-CPK41



MN-CPK42




MD-CPK42




MW-CPK42


C=4




MN-CPK43




MD-CPK43




MW-CPK43



MN-CPK44



MD-CPK44




MW-CPK44


NARROW

STANDARD


WIDE




MN-CPK61




MD-CPK61




MW-CPK61



MN-CPK62




MD-CPK62




MW-CPK62


C=6




MN-CPK63




MD-CPK63




MW-CPK63



MN-CPK64



MD-CPK64



MW-CPK64



NARROW


STANDARD

WIDE

STANDARD

WIDE

C=6



MN-CPC62




MD-CPC62



MW-CPC62



MN-CPC63



MD-CPC63



MW-CPC63



MN-CPC64




MD-CPC64




MW-CPC64


C=8




MD-CPK81




MW-CPK81




MD-CPK82




MW-CPK82




MD-CPK83



MW-CPK83



MD-CPK84



MW-CPK84

Cemented Restoration

CPK - Kit Contents




Tightening
force:
30Ncm




Burnout plastic cap
anti-rotation




Burnout plastic cap
free-rotation



Plastic healing cap



Transfer
plastic cap



Analog

CPK abutments are
also sold in kits

MK-NPK41 ● Ø4, H=1, C=4
MK-CPK41 ● Ø4.8, H=1, C=4
MK-WPK41 ● Ø5.5, H=1, C=4

MK-NPK42 ● Ø4, H=2, C=4
MK-CPK42 ● Ø4.8, H=2, C=4
MK-WPK42 ● Ø5.5, H=2, C=4

MK-NPK43 ● Ø4, H=3, C=4
MK-CPK43 ● Ø4.8, H=3, C=4
MK-WPK43 ● Ø5.5, H=3, C=4

MK-NPK44 ● Ø4, H=4, C=4
MK-CPK44 ● Ø4.8, H=4, C=4
MK-WPK44 ● Ø5.5, H=4, C=4

MK-NPK61 ● Ø4, H=1, C=6
MK-CPK61 ● Ø4.8, H=1, C=6
MK-WPK61 ● Ø5.5, H=1, C=6



MK-NPC62 ● Ø4, H=2, C=6
MK-CPC62 ● Ø4.8, H=2, C=6
MK-WPC62 ● Ø5.5, H=2, C=6

MK-NPK62 ● Ø4, H=2, C=6
MK-CPK62 ● Ø4.8, H=2, C=6
MK-WPK62 ● Ø5.5, H=2, C=6



MK-NPC63 ● Ø4, H=3, C=6
MK-CPC63 ● Ø4.8, H=3, C=6
MK-WPC63 ● Ø5.5, H=3, C=6

MK-NPK63 ● Ø4, H=3, C=6
MK-CPK63 ● Ø4.8, H=3, C=6
MK-WPK63 ● Ø5.5, H=3, C=6



MK-NPC64 ● Ø4, H=4, C=6
MK-CPC64 ● Ø4.8, H=4, C=6
MK-WPC64 ● Ø5.5, H=4, C=6

MK-NPK64 ● Ø4, H=4, C=6
MK-CPK64 ● Ø4.8, H=4, C=6
MK-WPK64 ● Ø5.5, H=4, C=6

MK-CPK81 ● Ø4.8, H=1, C=8
MK-WPK81 ● Ø5.5, H=1, C=8

MK-CPK82 ● Ø4.8, H=2, C=8
MK-WPK82 ● Ø5.5, H=2, C=8

MK-CPK83 ● Ø4.8, H=3, C=8
MK-WPK83 ● Ø5.5, H=3, C=8

MK-CPK84 ● Ø4.8, H=4, C=8
MK-WPK84 ● Ø5.5, H=4, C=8



Overdenture OT-Equator Kits

i Tightening
force:
30Ncm

NARROW



MK-NOE1




MK-NOE2



MK-NOE3



MK-NOE4



MK-NOE5

STANDARD



MK-SOE1




MK-SOE2



MK-SOE3



MK-SOE4



MK-SOE5

WIDE



MK-WOE1



MK-WOE3



MK-WOE5

Overdenture

OT-Equator Components



Tightening
force:
30Ncm

OT-Equator
handling tools



ET-IT002

OT-Equator impression
coping cap




OE-IC002

OT-Equator impression
transfer cap



OE-IC001

OT-Equator analog



OE-RSM60

Metal housing &
transfer caps



OE-MH001 OE-RCB01
OE-RCW01 OE-RCP01
OE-RCY01 OE-RCV01

OT-Equators
are sold in kits




Overdenture


Ball Attachment System

i Tightening force: 30Ncm

NARROW



STANDARD



WIDE



Overdenture

Ball Attachment System Components



Tightening
force:
30Ncm

Analog




MB-RS375

Gold / Plastic /
titanium ball cap



MB-GPS10 MB-TBS20
MB-SF200 MB-TBS10
MB-TTP10 MB-PR010

Metal housing &
Transfer caps



MB-DMH10
MB-PPC10
MB-PPR10

Plastic disc for ball
attachment



MB-DB235

Ball attachment
anchor key



MT-NB225

Ball attachment
anchor key



MT-RB225

Overdenture LOCKiT Attachment System

i Tightening
force:
30Ncm



NARROW



MN-LOC01



MN-LOC02



MN-LOC03



MN-LOC04



MN-LOC05

STANDARD



MD-LOC01



MD-LOC02



MD-LOC03



MD-LOC04



MD-LOC05

WIDE



MW-LOC01



MW-LOC03



MW-LOC05



All LOCKiT products shown in
this catalog are sold subject
to individual regulatory
certification only.

Overdenture LOCKiT System Components

i Tightening
force:
30Ncm

Abutment holder sleeve



MT-LOC11

Caps range




MM-BLA10 MM-CLE10
MM-BLU10 MM-RED10
MM-PIN10 MM-GRE10

Male processing package



MM-LOC02

Metal housing



MM-TIT10

Impression coping



MM-LOC11

LOCKiT analog



MM-LOA40
MM-LOA50

LOCKiT core tool




MT-LOC10

LOCKiT long ratchet
driver



MT-RSD24

LOCKiT short ratchet
driver




MT-RSD19

LOCKiT abutments
are also sold in kits



MK-NLC1	●	Ø3.75, H=1
MK-SLC1	●	Ø3.75, H=1
MK-WLC1	●	Ø4.5, H=1
MK-NLC2	●	Ø3.75, H=2
MK-SLC2	●	Ø3.75, H=2
MK-NLC3	●	Ø3.75, H=3
MK-SLC3	●	Ø3.75, H=3
MK-WLC3	●	Ø4.5, H=3
MK-NLC4	●	Ø3.75, H=4
MK-SLC4	●	Ø3.75, H=4
MK-NLC5	●	Ø3.75, H=5
MK-SLC5	●	Ø3.75, H=5
MK-WLC5	●	Ø4.5, H=5



Overdenture Locator Attachment System

i Tightening
force:
30Ncm


NARROW



STANDARD



WIDE



Overdenture Locator System Components

i Tightening
force:
30Ncm

Analogs



MM-LFA40
MM-LFA50

Metal housing &
transfer caps



MM-LMP02

Impression coping



MM-LIC10

Caps range:



MM-LLR10



MM-LER10



MM-LELE1



MM-LRM10



MM-LEL10



MM-LPP10



MM-LBP10

Tools:



MM-LCT10



MM-RSD15



MM-LSD15



MM-RSD20



MM-LSD21

Prosthetic Screws



Tightening
force:
20/30Ncm

Mini prosthetic screw




MD-S0200

Direct prosthetic screw




MN-S0160
MD-S0220

Medium prosthetic screw




MD-S0222

Long prosthetic screw



MD-S0224


Direct gold prosthetic screw



MD-G0220

NOTES





All rights reserved. No part of this publication may be reproduced, transcribed, stored in an electronic retrieval system, translated into any language or computer language, or be transmitted in any form whatsoever, without the prior written consent of the publisher. Warning: MIS products referred to in this document should be used by licensed dentists only.

